JRI ORTHOPAEDICS LTD

18 Churchill Way,
31A Business Park,
Chapelton,
Sheffield,
S35 2PY, UK

Instructions for Use

JRI Instrumentation – Sterile and Non-Sterile

JRI Instrumentation – Sterile & Non-Sterile

Contents

English.....	3
Español.....	12
Português	21
Norsk.....	30

Important Information

Further reprocessing information for JRI reusable instrumentation is available in IFU F4030, External Instrument Sterilisation Instruction. Please read instructions for use and corresponding operative technique, prior to use in a clinical setting. The Surgeon should be familiar with the appropriate operative technique. Current versions of operative techniques can be downloaded from the JRI Ltd website.

Caution: Contains Cobalt

Cobalt is classified as a carcinogenic, mutagenic or reprotoxic (CMR) substance and could result in an undesired biological reaction. Risks resulting from the presence of cobalt in this device have been duly assessed and minimised as far as possible.

Intended Purpose

JRI Sterile and Non-Sterile Instruments are intended to be used by orthopaedic surgeons to facilitate the implantation and subsequent explantation of compatible JRI implants in the way described in the respective operative technique.

Use of JRI instrumentation in accordance with the IFU and relevant operative technique is essential for successful implantation.

JRI Instrument Trays are intended to facilitate the containment, storage, transportation, cleaning, and sterilisation of reusable instruments.

JRI also offer general surgical instruments to facilitate surgical procedures which are not related to a specific JRI implant system or surgical technique.

Expected Clinical Benefits

JRI Sterile and Non-Sterile Instrumentation is not expected to provide a direct clinical benefit. The expected clinical benefit is derived from the associated JRI Implant.

Performance Characteristics

JRI Instrumentation is used to implant and subsequently explant compatible JRI implants. For information on the performance characteristics of JRI Implants, please refer to the IFU for the associated implant. All IFUs are available to download from the JRI Ltd website.

Treatment Population

For information on intended treatment populations, please refer to the IFU for the associated implant. All IFUs are available to download from the JRI Ltd website.

Description

Non-Sterile, Reusable Instrumentation

JRI Non-Sterile, Reusable Instrumentation includes a range of manual surgical instruments and instruments that may be used under power. These instruments are supplied non-sterile and are intended for sterilisation prior to each use by the end user. Validated reprocessing techniques are identified in JRI's Instrument Reprocessing IFU (F4030) which is available to download from the JRI website.

Sterile, Single-Use Instrumentation

JRI Sterile, Single-Use Instrumentation is sterilised using gamma irradiation and provided in a sterile state, as indicated on the instrument label. These instruments are intended to be used for a single patient in a single procedure only, and they may be used multiple times in this procedure. However, these instruments should NOT be resterilised or reused in multiple surgical procedures or on multiple patients. Reprocessing of these instruments may not result in adequate sterilisation and reuse may compromise the functionality and efficacy of the instrument.

If the sterile barrier packaging is broken, the instrument should NOT be used and cannot be re-sterilised for subsequent use.

Abbreviations and Definitions

12/14 = 12mm Small taper diameter, 14mm Large taper diameter

LS = Long Spigot. < = Angle AF = Across Flats

For Trial Stems: XXS = Extra Extra Small. XS = Extra Small. S = Small. M = Medium. L = Large. XL = Extra Large

For Trial Heads Neck lengths: S = Short. M = Medium. L = Long. XL = Extra Long. XL+4 = Extra Long plus 4mm. XL+8 = Extra Long plus 8mm

Warnings

- Non-sterile, reusable instruments **MUST** be cleaned and sterilised before use in accordance with the validated methods explained in JRI's Instrument Reprocessing IFU (F4030). **DO NOT** use devices which are suspected to be unclean or non-sterile at the point of use.
- Reuse of single use devices which have come into contact with blood, bone, tissue, or other body fluids may lead to injury. Potential risks associated with the reuse of single use devices include but are not limited to breakage, incorrect function, or failure of the device and transmission of infectious agents. These risks may result in patient injury or death.
- All instruments should be inspected for wear and damage prior to each use. In the event of visible wear or damage such as deformation, corrosion, or blunting of sharp/cutting instruments, the instrument **SHOULD NOT** be used as it may not perform as intended and risks injury of the patient and user.
- **DO NOT** subject instruments to excessive force or high loads as this may result in breakage. Where there is the requirement for the application of reasonable force, this will be explained in the operative technique.
- **DO NOT** modify or alter instruments unless explicitly instructed to do so in the operative technique or reprocessing IFU.
- **DO NOT** permanently implant any instrument. Ensure that instruments and parts of instruments **ARE NOT** left in the surgical site prior to closure, as this may result in injury to the patient. Instruments may not be detectable using commonly used imaging techniques such as radiography, fluoroscopy, and CT scanning.
- **DO NOT** use JRI Instruments in an MR environment. The instruments have not been tested for safety or compatibility in an MR environment. Use in such environment may result in heating or migration.
- Patients and Users with a hypersensitivity to the materials used in JRI Instruments, such as Cobalt or Nickel, may experience an allergic reaction.

Compatibility

JRI Instrumentation is intended for use with specific JRI implant systems. Instrumentation should only be used with other JRI compatible components. Instrument assemblies should only contain compatible components with corresponding connections. JRI Instruments and instrument components should not be used with other manufacturers' instruments unless specified in the operative technique.

Indications

For information on indications, please refer to the IFU for the associated implant. All IFUs are available to download from the JRI Ltd website.

Contra-Indications

For information on contra-indications, please refer to the IFU for the associated implant. All IFUs are available to download from the JRI Ltd website.

At a minimum, JRI instrumentation **MUST NOT** be used:

- Where there is active infection.

- By untrained, or inexperienced persons, or persons who are not clinically qualified to carry out the procedure for which the instruments are intended.
- For clinical procedures for which they are not intended.
- If they are visibly broken, or damaged intra-operatively, or fail to function as intended due to incorrect assembly.

Precautions

- Instruments must be carefully inspected before use to ensure they are free from scratches, dents, or other imperfections and damage that could result in breakage during use.
- In the event of breakage during use, the surgical site should be carefully inspected for instrument fragments to ensure that all fragments which may have entered the surgical site are removed prior to closure.
- Any decisions not to remove broken or fragmented instruments from the patient is at the discretion of the surgeon and must consider the associated risks. JRI accepts no liability for decisions not to remove broken or fragmented instruments.
- Non-sterile, reusable instruments must be cleaned, disinfected, sterilised, and inspected prior to each use.
- Any disinfection of JRI Instrumentation should be followed by sterilisation.
- Prolonged immersion in a disinfecting solution can damage surgical instruments therefore, they SHOULD NOT be soaked for a long period of time.
- Care should be taken during cleaning to ensure any holes, moving parts or other intricate areas of instrumentation are suitably cleaned and disinfected prior to sterilisation.
- Care should be taken when handling sharp instruments to reduce the risk of injury to the patient and user.
- Used instruments may be considered to be a biohazard and therefore, should be handled by competent personnel in accordance with local regulations and guidance.
- Care should be taken when handling and storing all instruments to prevent damage and breakage. As the manufacturer, JRI Orthopaedics Ltd can take no responsibility for damage, breakage or other adverse effects caused as a result of the failure of any person to follow these instructions or any other relevant applicable JRI instructions.
- Patients identified as risk patients with regard to Creutzfeldt-Jakob disease (CJD), related infections or other unusual transmissible diseases should be treated with single use instruments, which are disposed after surgery. Because of insufficient scientific evidence available today, a general efficient cleaning and decontamination procedure cannot be recommended. Current national recommendations should be considered.
- Lubrication oil is not required and should not be used on these instruments.

End Users

This instrumentation is to be used only under the control and supervision of an accredited Orthopaedic Surgeon or Physician whose responsibility it is to ensure that any user is qualified and trained in the use of these surgical instruments and the relevant surgical procedures. The medical team have a duty of care towards their patient which includes the correct use of this instrumentation. As a manufacturer of the product, JRI take no responsibility for any damage, breakage or adverse effects caused as a result of any failure in the medical team to discharge such duty.

Pre-Operative Preparation of Instruments

Sterile, single-use instruments are supplied sterile, and do not require cleaning, disinfection or sterilisation prior to use.

For non-sterile, reusable instrumentation, all packaging should be removed prior to decontamination and sterilisation, in accordance with the validated methods specified in the JRI Instrument Reprocessing IFU (F4030) and shall be maintained in a sterile state until used.

Pre-Operative

All instrumentation must be thoroughly inspected for possible damage prior to surgery.

The surgical staff should ensure that all reusable instrumentation has been adequately sterilised, using an approved, and validated, cleaning and decontamination process, and has been correctly reassembled prior to use, and is fit for purpose. For more complicated instrumentation, specific assembly Instructions are provided and should be consulted.

Clinicians should familiarise themselves with the instrumentation, prior to clinical use, particularly where that instrumentation requires intra-operatively assembly.

The surgeon should discuss all aspects of the surgery with the patient. Allergies and other reactions to instrument materials although rare should be considered and ruled out pre-operatively.

X-ray templates should be used to estimate implant sizes, placement and joint alignment.

Intra-Operative

Care should be taken not to cut through surgical gloves when handling any sharp-edged instruments and to consider the risk of infection if a cut appears.

JRI Instrumentation should only be used with the compatible JRI implant in accordance with the relevant operative technique. Operative Techniques are available to download from the JRI website. The surgeon is responsible for ensuring optimum implantation of the prosthetic device using JRI Instrumentation and should always refer to the relevant Operative Technique and Implant Instructions for Use.

Reasonable, but NOT excessive, force may be required in the use of this instrumentation. Where a specific force or torque is required, this is indicated in the operative technique.

Post-Operative

If the instrumentation is not fully intact or complete after the operation, then the surgeon should be adequately convinced that no parts are left *in vivo*. Fragments of instruments may be located by means of an image intensifier and/or X-ray radiography.

It is the responsibility of the surgical staff to ensure that all instrumentation is thoroughly cleaned and decontaminated as soon as possible after use.

Damaged instrumentation should be reported to your local JRI Sales Representative and returned to JRI Orthopaedic Ltd. for investigation.

Post-Operative Decontamination

Sterile, Single-Use instruments must not be resterilised, but be disposed of, in accordance with the corresponding hospital procedure for such components, following the completion of the surgical procedure.

Non-sterile, reusable instruments should be decontaminated immediately after use in accordance with F4030.

Any instrument that can be dismantled should have the various components cleaned separately, in accordance with the specific cleaning instructions provided with these instruments.

Attention should be paid to any cannulated instrument and any screw holes where debris may enter before being re-assembled. Recesses and hidden areas within an instrument should be inspected regularly to ensure that entrapped or other residual materials are completely removed.

If the instrument requires a complicated assembly technique, then refer to the additional assembly/dis-assembly instructions.

Where cleaned and decontaminated reusable instruments are to be returned to JRI, the appropriate decontamination certificate, provided in F4030, shall be completed by the Theatre Manager or their authorised Deputy

Adverse Effects

The following adverse effects may occur from the use of JRI Instrumentation:

Adverse Effects to Patient:

- Missing, Damaged, or worn instruments may cause surgical delays.
- Bone fractures may result from one-sided overloading or weakening of the bone substance. There may also be a risk of cardiovascular disorders, tissue reactions and haematoma.
- Use of blunt, damaged, or incorrectly positioned instruments may result in tissue injury or the need to remove additional bone.
- Improper reprocessing of instruments or use of sterile instruments where the sterility has been compromised may cause infection or toxicity.
- Sharp instruments may cause damage to vital organs.
- Sharp instruments or surgical trauma may cause neurological or vascular damage.

Adverse Effects to User:

- Sharp edges, burs, impaction, vibration or jamming of instruments may cause abrasions, cuts, contusions or other tissue injuries.

The surgeon is responsible for any complications that may result from erroneous indication, incorrect operative technique or inadequate aseptic precautions. Cleaning and decontamination instructions apply.

Further details on potential adverse effects can be found in the relevant implant IFU, including the potential adverse effects that are general to any surgery. All IFUs are available to download from the JRI website.

Storage and Handling

Sterile, Single Use Instrumentation should be stored and transported in their original protective packaging. Non-sterile, Reusable Instrumentation should be stored and transported in their original protective packaging or instrument tray. If an instrument tray is used, it should be wrapped using a double wrap method. All JRI Instrumentation should be kept out of sunlight in a clean and dry atmosphere. The acceptable storage period depends on the storage conditions, which must be determined by the Hospital.

Surgical instruments are sensitive to damage. Even small scratches can increase wear and the risk of corrosion. Instruments should always be handled with care.

Care should be taken when lifting complete instrument sets or transit cases.

Failure to follow these instructions for use may result in damage to or breakage of instruments and have potential adverse effects on patients and users. Misuse of instruments may reduce their useful life and increase the risk of injury to patients and users.

Device Disposal

Any disposal of JRI Instrumentation should follow hospital waste processing procedures.

Instrument Shelf-life

Non-sterile, reusable, instruments do not have a shelf-life. They are manufactured from materials that do not deteriorate over time, when stored unused under ambient conditions.

The shelf-life of sterile, single-use, instruments is indicated on the corresponding labelling for that instrument.

Functional Device Lifespan

JRI Instruments do not have an indefinite functional life. Many factors can affect the lifetime of JRI Instrumentation such as stresses related to surgical use, routine cleaning, and sterilisation processes. Visible signs of wear to the instruments could indicate that the instrument may not function as expected. Therefore, all instrumentation is required to be inspected for signs of wear, such as scratches, dents, and other damage, before each use, and, where appropriate, reported to JRI Orthopaedics Ltd, for reprocessing, or replacement. The cutting edges of sharp instruments should be maintained and replaced at the surgeons' discretion.

MRI Safety Information

JRI Instruments are not intended for use in an MR environment. They have not been tested for safety or compatibility in an MR environment. The instruments have not been tested for heating or migration in an MR environment.

Further Information

For product specific operative technique training or further information, please contact your JRI Orthopaedics Ltd. Sales Representative or JRI Orthopaedics Ltd. directly.

Adverse Event Reporting

Any serious incident that has occurred in relation to the devices should be reported to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Incidents can be reported to the manufacturer using the following email: regulatory@jri-ltd.co.uk

JRI ORTHOPAEDICS LTD

18 Churchill Way,
35A Business Park,
Chapelton,
Sheffield,
S35 2PY, UK
Tel: +44(0)114 345 0000
Fax: +44(0)114 345 0004
www.jri-ltd.co.uk

EC	REP
----	-----

CMC Drugs & Medical Devices S.L.
C/Horacio Lengo No. 18,
CP29006,
Malaga, Spain
Tel: +3495124054
Fax: +34592330100
www.cmcmedicaldevices.com

Symbol Glossary

Applicable Symbols	Description
	eIFU Indicator. Electronic IFU is available on the manufacturer's website. Consult the instructions for use.
	Keep away from sunlight.
	Keep dry.
	Do not use if packaging is damaged and consult instructions for use.
	Single-use device. Do not re-use.
	Caution – check for specific warnings or precautions.
	Non-sterile.
	Sterilised using irradiation.
	Double sterile barrier system.
	Single sterile barrier system with protective packaging inside.
	Batch number, batch code or lot number.
	Catalogue number.

	<p>Unique device identifier.</p>
	<p>Manufacturer.</p>
	<p>Date of manufacture.</p>
	<p>Use-by date.</p>
	<p>Do not resterilise.</p>
	<p>Medical device.</p>
	<p>Authorised representative in the European Community.</p>
	<p>The device meets the essential requirement of the European Medical Device Regulations.</p>
	<p>Contains hazardous substances.</p>

Revision History

Issue	Date	Change
12	07 Mar 2025	<p>Updated, as follows, to meet MDR requirements: Removed the materials section. Added the following sections:</p> <ul style="list-style-type: none"> • Caution: Contains Cobalt • Expected Clinical Benefit • Precautions • Warnings • End Users • Disposal • MRI Safety Information • Adverse Event Reporting • Symbol Glossary • Revision History <p>Updated the following sections:</p> <ul style="list-style-type: none"> • Intended Purpose • Treatment Population • Description • Pre-Operative Preparation of Instruments • Pre-Operative • Intra-Operative • Post-Operative Decontamination • Adverse Effects • Storage And Handling • Functional Device Lifespan <p>Updated branding.</p>

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